

■ Company Profile

Name of Company : MITAKA Co., Ltd.

President : Eigo Tajima

Head Office : 1765-1, Menuuma, Kumagaya-shi,
Saitama, 360-0201, Japan

Date of Establishment : October, 1991

Capital : 10,000,000JPY

Factory #1 : 1765-1, Menuuma, Kumagaya-shi,
Saitama, 360-0201, Japan

TEL : +81-48-588-5552

FAX : +81-48-589-0040

Factory #2 : 638-8, Yatogo, Kumagaya-shi, 360-0203,
Japan

Major Banks : Gunma Bank

■ Company history

1991 Oct. Establishment of MITAKA Co., Ltd.

1994 July. Change of registered head office and
factory address to Yatogo, Kumagaya-shi

1998 Oct. Construct new factory at Yatogo,
Kumagaya-shi

2001 Jan. Establishment of Dong Guan Chang Ping
Kyoritsu Metal Factory

2003 March. Establishment of M.T.K. Technology
Company Limited in Hong Kong

2010 Feb. Change of registered head office and
factory address to Menuuma, Kumagaya-shi

2011 March Awarded ISO9001

2012 Jan. Capital increase to 10,000,000JPY

■ ISO Certification

ISO9001 : 2015


■ Process technology

● Coiling Micro coils

Provide precious metal coils with high precision
and excellent response to tip coil and radiopaque
marker for Guide Wire

Material : Pt Alloy Au SUS(Medical)、Music
Wire(Industrial) etc

Wire Diameter : 10 μ m~



● Coiling Long coils

Provide high accuracy and high operability long
coils for protective sheaths of endoscopic treatment
tools and various kind of catheters

Material : SUS(Medical)


Wire Diameter : 50 μ m~ Round wire, Flat wire



● Straightened and Cut Wire

Provide straightened wires excellent in
straightness and torque transmission performance to
medical Guide Wire

Material : SUS(Medical)、Copper Alloy(Industrial)
Wire Diameter : 20 μ m~



● Forming Grinding Polishing

Form the wire in various shapes, grind a part of the
wire into a taper shape, and polish the terminal to a
hemisphere

Material : SUS(Medical)、Tungsten(Industrial)


Wire Diameter : 0.1 mm~



■ Business items


● Medical coils

Platinum coil, Precious metal coil, Long coil.



● Medical parts

Stylet, Mandrel, Core wire.



● Others

Polygon coil, PTFE coating, etc.

■ Access

● By Train and Bus


From Kumagaya Station by JR Takasaki Line, get on Asahi Motor Corporation Bus heading to Ota Station, get off at Menuma-shotenmae and walk 10min.

● By Car

Kan-Etsu Expressway exit at Higashimatsuyama Interchange, drive for 26Km

Tohoku Expressway exit at Hanyu Interchange, drive for 39Km

● Map



■ Homepage

<http://www.t-mtk.co.jp/original10.html>

■ Youtube

<https://youtube.com/channel/UCUNU7IvB-9SBDqPtau9Lz3w>

■ Contact us

e-mail : info@t-mtk.co.jp

Tel : +81-48-588-5552

Contribute to metal wire processing technology of
the latest medical industry

